





**LENUS MEDICARE & RESEARCH Pvt.Ltd.**


[www.lenusmedicare.com](http://www.lenusmedicare.com)



## MASTER STENT PUSHER

A biliary master stent pusher is a precise medical tool that helps place biliary stents correctly during endoscopic or percutaneous procedures. It ensures smooth delivery, control, and positioning of the stent for effective biliary drainage. This device is perfect for use in endoscopic retrograde cholangiopancreatography (ERCP) or percutaneous biliary treatments, providing excellent performance and dependability in both diagnosis and treatment. It is compatible with a wide range of biliary stents, ensuring versatility across various clinical applications.

This medical device is used to drain obstructed biliary ducts. The characteristics of the Biliary Plastic Stent Introducer Set are radiopaque material for 10F guiding catheter and a smooth surface that helps in gliding smoothly over the guidewire. The device is sterilized, ready-to-use, and available in multiple sizes to suit procedural needs.



# PRODUCT SPECIFICATION

## SPECIFICATION DETAILS

| Length (cm) | Diameter (Fr) | Guidewire(in) | Minimum working Channel (mm) | Wall thickness (catheter) (mm) | Marker length (mm) |
|-------------|---------------|---------------|------------------------------|--------------------------------|--------------------|
| 180         | 10            | 0.035"        | 3.7                          | 1.14                           | 3                  |

## KEY FEATURES

**Accurate Stent Placement:** Provides precise control for optimal positioning of biliary stents.

**Smooth Delivery System:** Minimizes friction for effortless stent deployment.

**Ergonomic Design:** Ensures comfortable grip and easy maneuverability.

**Biocompatible Materials:** Made from medical-grade, non-toxic components for patient safety.

**Sterile & Ready to Use:** Pre-sterilized for immediate clinical use.

**Compatible with Various Stents:** Works with multiple biliary stent types and sizes.

**Radiopaque Markers :** Enhances visibility under fluoroscopy for accurate placement.

## PRODUCT FEATURES & GUIDE

### PROCEDURE GUIDE

Insert the stent pusher carefully through the biopsy channel.

It passed down your throat and through the stomach to the small intestine.

The camera shows the inside of your body on a screen.

The stent is put inside the bile duct using a small wire.

When the stent is in the right place the wire is removed.

Withdraw the catheter from the endoscope.